

Beta Glucosidase, Leukocytes

Order Name: **BETA GLUCS**

Test Number: 3661600

Revision Date: 12/09/2021

| TEST NAME | METHODOLOGY | LOINC CODE |
|------------------------------|-------------|------------|
| Beta Glucosidase, Leukocytes | | 32540-7 |

SPECIMEN REQUIREMENTS

| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
|-----------|-----------------------|---------------|----------------------------------|-----------------------|
| Preferred | 6.0mL (5.0mL) | Whole Blood | ACD Solution A or B (Yellow Top) | Refrigerated |

Instructions

For optimal isolation of leukocytes, it is recommended the specimen arrive refrigerated within 6 days of collection to be stabilized. Collect specimen **Monday through Thursday only and not the day before a holiday.**

Specimen should be collected and packaged as close to shipping time as possible.

Unacceptable: Gross hemolysis

Specimen Stability: Room temperature: 6 days, Refrigerated: 6 days, Frozen: unacceptable

Submit specimen in original collection container.

Please submit the following forms with patient specimen:

[Consent for Genetic Testing Form](#)

[Mayo Patient Information Form](#)

GENERAL INFORMATION

| | |
|------------------|--|
| Testing Schedule | Monday |
| Expected TAT | 8-15 days |
| Notes | Reference Lab: Mayo Mayo Test Code: GBAW Click Here to view information on the Mayo website. |
| CPT Code(s) | 82963 |
| Lab Section | Reference Lab |