

Herpes Select Panel (HSV 1/2 IgG and IgM)

Order Name: **Herpes Panel**

Test Number: 5563986

Revision Date: 04/20/2020

| TEST NAME | METHODOLOGY | LOINC CODE |
|---------------------------------------|-------------------------------|------------|
| HSV Type 1 IgG | Chemiluminescence Assay | 5206-8 |
| HSV Type 1 IgG Interpretation | | 51916-5 |
| HSV Type 2 IgG | Chemiluminescence Assay | 5209-2 |
| HSV Type 2 IgG Interpretation | | 43180-9 |
| Herpes Simplex Virus Type 1/2 IgM IFA | Indirect Fluorescent Antibody | |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|------------------------------------------------------------------------------------------|---------------|--------------------|------------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 3 mL | Serum | Clot Activator SST | Refrigerated or Frozen |
| Instructions | Specimen Stability: Room Temperature: Unacceptable, Refrigerated: 7 days, Frozen 30 Days | | | |

| GENERAL INFORMATION | |
|---------------------|----------------------------------------------------------------------------|
| Testing Schedule | Varies |
| Expected TAT | 3 - 5 days |
| Notes | HSV 1/2 IgG are performed Mon - Sat and HSV 1/2 IgM are performed Tues/Fri |
| CPT Code(s) | 86695x2, 86696x2 |
| Lab Section | Immunology - Serology |