Labcorp Oklahoma, Inc. Test Directory

Biotinidase, Serum (with Paired Normal Control)

Order Name:BiotinidaseTest Number:5613879Revision Date:12/12/2022

TEST NAME			METHODOLOGY	LOINC CODE
Biotinidase, Serum (with Paired Normal Control)			Spectrophotometry	
SPECIMEN REQUIR	REMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (1 mL)	Serum	Clot Activator (Red Top, No-Gel)	Frozen
Instructions	Specimen Type: Red-top tube Specimen Storage: FREEZE Specimen Collection: Separate serum from cells within 30 minutes of venipuncture. Transfer specimen to plastic transport tube before freezing. To avoid delays in turnaround time when requesting multiple tests on frozen samples, please submit separate frozen specimens for each test requested. Special Instructions: ***NOTE - DO NOT USE THIS TEST CODE FOR BIOTIN/VIT B7*** TESTING Do not modify test inquiry names without permission of CMBP.			

GENERAL INFORMATION		
Expected TAT	4-10 days	
Performing Labcorp Test Code	402362	
Notes	Labcorp Test Code: 402362	
CPT Code(s)	82261	
Lab Section	Reference Lab	