

AH50 (Complement Activity, Alternative Pathway)

Order Name: **AH50**
Test Number: 5002275
Revision Date: 04/14/2023

| TEST NAME | METHODOLOGY | LOINC CODE |
|---|-------------|------------|
| AH50 (Complement Activity, Alternative Pathway) | Hemolysis | 13087-2 |

SPECIMEN REQUIREMENTS

| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
|--------------|--|---------------|----------------------------------|-----------------------|
| Preferred | 1 mL | Serum | Clot Activator (Red Top, No-Gel) | Frozen |
| Instructions | Specimen Type: Red-top tube Specimen Storage: FREEZE (SHIP ON DRY ICE) Do Not Share this specimen with another test. Specimen Collection: Transfer specimen to plastic transport tube before freezing. Send to lab on dry ice. Special Instructions: Testing referred to National Jewish Health NJHCO#AM50 Specimen Stability: Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available | | | |

GENERAL INFORMATION

| | |
|------------------------------|---------------------------|
| Expected TAT | 27 - 36 days |
| Performing Labcorp Test Code | 831782 |
| Notes | Labcorp Test Code: 831782 |
| CPT Code(s) | 86161 |
| Lab Section | Reference Lab |