

Candida Vaginitis - Swab - TMA

Order Name: **Cand Vag Swab**

Test Number: 6987000

Revision Date: 04/06/2020

| TEST NAME | METHODOLOGY | LOINC CODE |
|--------------------------------|--------------------------------------|------------|
| Candida Species - Swab | Transcription-Mediated Amplification | 47000-5 |
| Candida Glabrata - Swab | Transcription-Mediated Amplification | n/a |

SPECIMEN REQUIREMENTS

| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|----------------------------------|-------------------------------|
| Preferred | Swab | Swab | APTIMA Multitest Swab kit | Ambient / Refrigerated |
| Instructions | Collect Swab in APTIMA Multitest Swab collection kit. Place swab in Swab Specimen Transport Tube, break shaft off at scoreline then recap tube. Stability: Room Temperature 30 Days, Refrigerated 30 Days, Frozen: N/A | | | |

GENERAL INFORMATION

| | |
|-------------------------|-----------------------|
| Testing Schedule | Mon-Fri |
| Expected TAT | 1-2 days after set-up |
| CPT Code(s) | 87481 |
| Lab Section | Molecular Diagnostic |