

Anemia General Panel

Order Name: **Anemia Gen Pnl**

Test Number: 2930825

Revision Date: 06/17/2021

| TEST NAME            | METHODOLOGY | LOINC CODE |
|----------------------|-------------|------------|
| Anemia General Panel |             |            |

| SPECIMEN REQUIREMENTS |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |                  |                                            |                       |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|--------------------------------------------|-----------------------|
| Specimen              | Specimen Volume (min)                                                                                                                                                                                                                                                                                                                                                                                                                                                           | Specimen Type    | Specimen Container                         | Transport Environment |
| Preferred             | See Instructions                                                                                                                                                                                                                                                                                                                                                                                                                                                                | See Instructions | EDTA (lavender top) and Clot Activator SST | Refrigerated          |
| Instructions          | <p><b>Collect Both:</b> One 5mL(3mL) EDTA (Lavender) and One 5mL Clot Activator SST.</p> <p>For best results: <b>Room temperature specimens should be tested within 24hrs, otherwise send Refrigerated.</b></p> <p>Refrigerated specimens can be tested up to 48hrs. Specimens received after 24hrs will not receive a manual differential.</p> <p>Slides must be made within 24hr of collection.</p> <p><b>Specimens received greater than 48hrs old will be canceled.</b></p> |                  |                                            |                       |

| GENERAL INFORMATION |                            |
|---------------------|----------------------------|
| Testing Schedule    | Daily                      |
| Expected TAT        | 1 Day                      |
| CPT Code(s)         | 85025, 82607, 83550, 83540 |
| Lab Section         | Hematology                 |