

## VMA, Random Urine

Order Name: **VMA U Rand**

Test Number: 5194550

Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
VMA, Random Urine		

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
<b>Instructions</b>	<b>Specimen Type:</b> Plastic urine container <b>Specimen Storage:</b> Vanillylmandelic acid is STABLE at room temperature unpreserved for seven days. <b>Specimen Collection:</b> If specimen is to be stored for longer than seven days before analysis, use random urine transport tube. <b>Special Instructions:</b> Use for patients where 24-hour urine collection is difficult (eg, pediatric patients). Twenty-four hour urine collection is preferable (see Vanillylmandelic Acid (VMA), 24-Hour Urine [004143]). <b>Specimen Stability:</b> Ambient: 7 days, Refrigerated : 14 days, Frozen: 14 days			

### GENERAL INFORMATION

<b>Performing Labcorp Test Code</b>	123208
<b>Notes</b>	Labcorp Test Code: 123208
<b>CPT Code(s)</b>	82570, 84585
<b>Lab Section</b>	Reference Lab