

VMA, Random Urine

Order Name: **VMA U Rand**
Test Number: 5194550
Revision Date: 12/12/2022

| TEST NAME | METHODOLOGY | LOINC CODE |
|-------------------|-------------|------------|
| VMA, Random Urine | | |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|---|---------------|--------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Instructions | <p>Specimen Type: Plastic urine container</p> <p>Specimen Storage: Vanillylmandelic acid is STABLE at room temperature unpreserved for seven days.</p> <p>Specimen Collection: If specimen is to be stored for longer than seven days before analysis, use random urine transport tube.</p> <p>Special Instructions: Use for patients where 24-hour urine collection is difficult (eg, pediatric patients). Twenty-four hour urine collection is preferable (see Vanillylmandelic Acid (VMA), 24-Hour Urine [004143]).</p> <p>Specimen Stability: Ambient: 7 days, Refrigerated : 14 days, Frozen: 14 days</p> | | | |

| GENERAL INFORMATION | |
|------------------------------|---------------------------|
| Performing Labcorp Test Code | 123208 |
| Notes | Labcorp Test Code: 123208 |
| CPT Code(s) | 82570, 84585 |
| Lab Section | Reference Lab |