

Flunitrazepam (Rohypnol) Scr.U

Order Name: **FNZ Scr Ur**

Test Number: 5194624

Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Flunitrazepam (Rohypnol) Scr.U	Immunoassay (IA)	

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	10 mL (3 mL)	Urine, Random	Sterile Urine container	Room Temperature
Instructions	<p>Notes: 3 mL (Note: This volume Does NOT allow for repeat testing). Specimen Type: Plastic urine container without preservative Specimen Storage: SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Special Instructions: Testing referred to MedTox Laboratories Inc MX#436</p>			

GENERAL INFORMATION

Expected TAT	3-10 Days
Performing Labcorp Test Code	808417
Notes	Labcorp Test Code: 808417
CPT Code(s)	80307
Lab Section	Reference Lab