

Venlafaxine/Metabolite Level

Order Name: **VENLAFAXIN**

Test Number: 4503150

Revision Date: 12/30/2022

| TEST NAME | METHODOLOGY | LOINC CODE |
|------------------------------|--|------------|
| Venlafaxine/Metabolite Level | High-pressure liquid chromatography (HPLC) with fluorometric detection | |

SPECIMEN REQUIREMENTS

| Specimen | Specimen Volume (mL) | Specimen Type | Specimen Container | Transport Environment |
|--------------|--|---------------|---|-----------------------|
| Preferred | 3 mL (0.6 mL) | Plasma | Lithium Heparin (Dark Green Top / No-Gel) | Refrigerated |
| Alternate 1 | 3 mL (0.6 mL) | Serum | Clot Activator (Red Top, No-Gel) | Refrigerated |
| Instructions | <p>Specimen Type: Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended.</p> <p>Specimen Storage: SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.</p> <p>Specimen Collection: Serum or plasma should be separated from cells within two hours of venipuncture. Submit serum or plasma in a plastic transport tube.</p> <p>Special Instructions: Testing referred to MEDTOX Laboratories Inc TC 179</p> | | | |

GENERAL INFORMATION

| | |
|------------------------------|---------------------------|
| Expected TAT | 7 - 10 days |
| Performing Labcorp Test Code | 808514 |
| Notes | Labcorp Test Code: 808514 |
| CPT Code(s) | 80299 |
| Lab Section | Reference Lab |