Labcorp Oklahoma, Inc. Test Directory

Ethanol Biomarkers, MS, Ur, MAT

Order Name: Etoh Biomark Ur

Test Number: 5195577 Revision Date: 01/09/2024

TEST NAME			METHODOLOGY	LOINC CODE
Ethanol Biomarkers,MS, Ur, MAT		Mass Spectrometry (MS)		
SPECIMEN REQU	IREMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	30 mL (15 mL)	Urine, Random	Sterile Urine Container Preservative Free	Room Temperature
Instructions	Specimen: 30mL(15mL) Random Urine Preservative free Urine Transport tube Room Temperature or Refrigerated. Cause for Rejection: Insufficient volume, no ID on container, urine from preservative tube. SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.			

GENERAL INFORMATION		
Expected TAT	5-8 days	
Performing Labcorp Test Code	701737	
CPT Code(s)	80321	
Lab Section	Reference Lab	

Stability: Room Temperature 5 days, Refrigerated 14 days, Frozen 14 days.

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025