

ToxAssure Flex 23, Urine

Order Name: **TA Flex 23 UR**
 Test Number: 5195609
 Revision Date: 03/31/2025

TEST NAME	METHODOLOGY	LOINC CODE
ToxAssure Flex 23, Urine	Immunoassay (IA)	

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	30 mL (3 mL)	Urine, Random	Sterile Urine Container Preservative Free	Room Temperature
Instructions	Specimen: 30mL(3mL) Random Urine Preservative free Urine Transport tube Room Temperature or Refrigerated. Cause for Rejection: Insufficient volume, no ID on container, urine from preservative tube. SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Stability: Room Temperature 14 days, Refrigerated 14 days, Frozen 6 weeks.			

GENERAL INFORMATION

Expected TAT	5-8 days
Performing Labcorp Test Code	912335
CPT Code(s)	80362, 80348, 80368, 80347, 80357, 80354, 80307
Lab Section	Reference Lab