## Labcorp Oklahoma, Inc. Test Directory

## **Phenobarbital Level**

Order Name: PHENOBARB

Test Number: 4003300 Revision Date: 05/28/2021

| TEST NAME             |  |               | METHODOLOGY                           | LOINC CODE            |  |
|-----------------------|--|---------------|---------------------------------------|-----------------------|--|
| Phenobarbital Level   |  |               | Enzyme Immunoassay                    | 3948-7                |  |
| SPECIMEN REQUIREMENTS |  |               |                                       |                       |  |
| Specimen              | Specimen Volume (min)  | Specimen Type | Specimen Container                    | Transport Environment |  |
| Preferred             | 1 mL (0.5)   | Plasma        | Lithium Heparin PST (Light Green Top) | Refrigerated          |  |
| Alternate 1           | 1 mL (0.5)   | Serum         | Clot Activator SST                    | Refrigerated          |  |
| Instructions          | Draw specimen same time each day. Specimen stability: Ambient 7 days, Refrigerated 6 months, Frozen 6 months |               |                                       |                       |  |

| GENERAL INFORMATION |                     |  |
|---------------------|---------------------|--|
| Testing Schedule    | Daily               |  |
| Expected TAT        | 1-2 days            |  |
| Notes               | For toxicology use. |  |
| CPT Code(s)         | 80184               |  |
| Lab Section         | Chemistry           |  |
|                     |                     |  |

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025