## Flow Cytometry on Body Fluid

Order Name: FLOW FLUID
Test Number: 8090085
Revision Date: 07/06/2022

TEST NAME			METHODOLOGY	LOINC CODE	
Flow Cytometry on Body Fluid		Flow cytometry			
SPECIMEN REQUIR	EMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	5 mL (3 mL)	Body Fluid	Sterile Screwtop Container	Refrigerated	
Instructions	Aliquot CSF, Thoracentesis/Pleural fluid or Paracentesis/Ascites fluid in a sterile screwtop container.  Keep Specimen Refrigerated. Indicate type of fluid on specimen tube and requisition.  Please deliver to performing laboratory department (flow cytometry) ASAP.  Specimen stability: 48hrs refrigerated				
OENEDAL INFORM	ATION				
GENERAL INFORMATION					

GENERAL INFORMATION		
Testing Schedule	Tue-Sat	
Expected TAT	2 Days	
CPT Code(s)	test dependant	
Lab Section	Flow Cytometry	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2024