Labcorp Oklahoma, Inc. Test Directory

Beta Glucosidase, Leukocytes

Mayo Patient Information Form

Order Name: **BETA GLUCS**

Test Number: 3661600 Revision Date: 12/09/2021

TEST NAME			METHODOLOGY	LOINC CODE	
Beta Glucosidase, Leukocytes				32540-7	
SPECIMEN REQU	JIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	6.0mL (5.0mL)	Whole Blood	ACD Solution A or B (Yellow Top)	Refrigerated	
Instructions	For optimal isolatoin of leukocytes, it is recommended the specimen arrive refrigerated within 6 days of collection to be stabilized. Collect specimen Monday through Thursday only and not the day before a holiday.				
	Specimen should be collected and packaged as close to shipping time as possible. Unacceptable: Gross hemolysis				
	Specimen Stability: Room temperature: 6 days, Refrigerated: 6 days, Frozen: unacceptable				
	Submit specimen in original collection container.				
	Please submit the following forms with patient specimen:				
	Consent for Genetic Testing Form				

GENERAL INFORMATION		
Testing Schedule	Monday	
Expected TAT	8-15 days	
Notes	Reference Lab: Mayo Mayo Test Code: GBAW Click Here to view information on the Mayo website.	
CPT Code(s)	82963	
Lab Section	Reference Lab	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025