

Gonadotropin releasing hormone (aka: LHRH)

Order Name: **Gonadotrop RH**

Test Number: 6905995

Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Gonadotropin releasing hormone (aka: LHRH)	High-pressure liquid chromatography (HPLC) tandem mass spectrometry.	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	3 mL (0.5 mL)	Plasma	EDTA (Lavender Top)	Frozen
Instructions	<p>Notes: 0.5 mL (Note: This volume Does NOT allow for repeat testing)</p> <p>Specimen Type: Trasylol(R) collection kit (Obtain no 33328 through LabCorp branch locations)</p> <p>Specimen Storage: FREEZE</p> <p>Specimen Collection: Using a chilled lavender-top (EDTA) tube, collect a whole blood specimen, add 0.25 mL of Trasylol(R) to the tube and mix well. Centrifuge the sample and remove the plasma. FREEZE immediately. To avoid delays in turnaround time when requesting multiple tests on frozen samples, PLEASE SUBMIT SEPARATE FROZEN SPECIMENS FOR EACH TEST REQUESTED.</p> <p>Special Instructions: Testing is performed at Esoterix Endocrinology Laboratory: ES#501836 NOTE: Add aprotinin (Trasylol(R)) to all specimens to allow for analyte stability. Separate plasma immediately, store, and ship frozen.</p> <p>Specimen Stability: Ambient: Not Available, Refrigerated : 48 hours, Frozen: 291 days</p>			

GENERAL INFORMATION	
Expected TAT	8 - 12 days
Performing Labcorp Test Code	502421
Notes	Labcorp Test Code: 502421
CPT Code(s)	83520
Lab Section	Reference Lab