

Western Waterhemp IgE

Order Name: **West Waterhemp**

Test Number: 6906843

Revision Date: 04/14/2020

TEST NAME	METHODOLOGY	LOINC CODE
W Waterhemp Count	ImmunoCAP	
W Waterhemp Class	ImmunoCAP	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (1 mL)	Serum	Clot Activator SST	Refrigerated
Instructions	Specimen: Collect 1 -2 ml of whole blood into a SST tube; centrifuge and transfer 0.5 ml of serum into a transfer tube Stability: Ambient - 4 weeks; Refrigerated - 4 weeks; Frozen - 4 weeks			

GENERAL INFORMATION	
Testing Schedule	Mon - Fri
Expected TAT	1 - 2 days
Notes	Reference lab: Viracor Test code: 72710S Click Here to view information on the Viracor website
CPT Code(s)	86003
Lab Section	Reference Lab