

Propoxyphene, MS, Ur, MAT

Order Name: **Propoxyphene UR**

Test Number: 5195585

Revision Date: 01/09/2024

| TEST NAME | METHODOLOGY | LOINC CODE |
|---------------------------|------------------------|------------|
| Propoxyphene, MS, Ur, MAT | Mass Spectrometry (MS) | |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|---|---------------|--|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 30 mL (15 mL) | Urine, Random | Sterile Urine Container Preservative Free | Room Temperature |
| Instructions | Specimen: 30mL(15mL) Random Urine Preservative free Urine Transport tube Room Temperature or Refrigerated. Cause for Rejection: Insufficient volume, no ID on container, urine from preservative tube. SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Stability: Room Temperature 5 days, Refrigerated 14 days, Frozen 14 days. | | | |

| GENERAL INFORMATION | |
|------------------------------|---------------|
| Expected TAT | 5-8 days |
| Performing Labcorp Test Code | 701746 |
| CPT Code(s) | 80367 |
| Lab Section | Reference Lab |