

## Pneumococcal Ab Pnl PCV20, 20 Serotypes

Order Name: **Pneumo 20**

Test Number: 5197339

Revision Date: 03/16/2026

TEST NAME	METHODOLOGY	LOINC CODE
Pneumococcal Ab Pnl PCV20, 20 Serotypes	Multiplex fluorescent immunoassay	See Below

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5 mL)	Serum	Clot Activator (SST or Red No-Gel)	Refrigerated

**Instructions**

**Specimen:** 1mL(0.5mL) Serum from SST clot tube  
Blood should be allowed to clot 30-60 minutes, then poured off into a plastic transport tube

**Storage Instructions:** Refrigerated

**Stability Requirements:** Ambient 6 weeks, Refrigerated 6 weeks, Frozen 6 weeks

**Causes for Rejection:** Insufficient sample volume.

Labcorp Oklahoma, Inc.  
Test Directory

**GENERAL INFORMATION**

<b>Expected TAT</b>	7-10 Days
<b>Clinical Use</b>	Pneumococcal antibody panel provides quantitative measurement of serotype-specific IgG antibodies to assess immune response to pneumococcal vaccination, particularly in immunocompromised patients.
<b>Performing Labcorp Test Code</b>	834975

**Notes**

RESULT CODE	RESULT NAME	RESULT LONG NAME	LOINC CODE
5575653	Pneumo 1 IgG	Pneumococcal Serotype 1 IgG	85955-3
5575659	Pneumo 3 IgG	Pneumococcal Serotype 3 IgG	86081-7
5197340	Pneumo 6A IgG	Pneumococcal Serotype 6A IgG	n/a
5575661	Pneumo 4 IgG	Pneumococcal Serotype 4* IgG	86108-8
5575663	Pneumo 5 IgG	Pneumococcal Serotype 5 IgG	86129-4
5575689	Pneumo 6B IgG	Pneumococcal Serotype 6B* IgG	40905-2
5575697	Pneumo 7F IgG	Pneumococcal Serotype 7F IgG	40911-0
5575667	Pneumo 8 IgG	Pneumococcal Serotype 8 IgG	86148-4
5575707	Pneumo 9V IgG	Pneumococcal Serotype 9V*, IgG	40926-8
5575691	Pneumo 10A IgG	Pneumococcal Serotype 10A IgG	86099-9
5575693	Pneumo 11A IgG	Pneumococcal Serotype 11A IgG	86123-7
5575671	Pneumo 12F IgG	Pneumococcal Serotype 12F IgG	85974-4
5575673	Pneumo 14 IgG	Pneumococcal Serotype 14* IgG	85992-6
5575699	Pneumo 15B IgG	Pneumococcal Serotype 15B IgG	47356-1
5575701	Pneumo 18C IgG	Pneumococcal Serotype 18C* IgG	40913-6
5575703	Pneumo 19A IgG	Pneumococcal Serotype 19A IgG	40915-1
5575679	Pneumo 19F IgG	Pneumococcal Serotype 19F* IgG	86021-3
5575683	Pneumo 22F IgG	Pneumococcal Serotype 22F IgG	86053-6
5575687	Pneumo 23F IgG	Pneumococcal Serotype 23F* IgG	86061-9
5575709	Pneumo 33F IgG	Pneumococcal Serotype 33F IgG	47349-6

**CPT Code(s)** 86581

**Lab Section** Reference Lab