

## Pneumococcal Ab Pnl PCV21, 21 Serotypes

Order Name: **Pneumo 21**

Test Number: 5197337

Revision Date: 03/16/2026

TEST NAME	METHODOLOGY	LOINC CODE
Pneumococcal Ab Pnl PCV21, 21 Serotypes	Multiplex fluorescent immunoassay	See Below

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5 mL)	Serum	Clot Activator (SST or Red No-Gel)	Refrigerated

**Instructions**

**Specimen:** 1mL(0.5mL) Serum from SST clot tube  
Blood should be allowed to clot 30-60 minutes, then poured off into a plastic transport tube

**Storage Instructions:** Refrigerated

**Stability Requirements:** Ambient 6 weeks, Refrigerated 6 weeks, Frozen 6 weeks

**Causes for Rejection:** Insufficient sample volume.

Labcorp Oklahoma, Inc.  
Test Directory

**GENERAL INFORMATION**

<b>Expected TAT</b>	7-10 Days
<b>Clinical Use</b>	Pneumococcal antibody panel provides quantitative measurement of serotype-specific IgG antibodies to assess immune response to pneumococcal vaccination, particularly in immunocompromised patients.
<b>Performing Labcorp Test Code</b>	835119

**Notes**

RESULT CODE	RESULT NAME	RESULT LONG NAME	LOINC CODE
5197340	Pneumo 6A IgG	Pneumococcal Serotype 6A IgG	n/a
5575659	Pneumo 3 IgG	Pneumococcal Serotype 3 IgG	86081-7
5575697	Pneumo 7F IgG	Pneumococcal Serotype 7F IgG	40911-0
5575667	Pneumo 8 IgG	Pneumococcal Serotype 8 IgG	86148-4
5575669	Pneumo 9N IgG	Pneumococcal Serotype 9N IgG	86166-6
5575691	Pneumo 10A IgG	Pneumococcal Serotype 10A IgG	86099-9
5575693	Pneumo 11A IgG	Pneumococcal Serotype 11A IgG	86123-7
5575671	Pneumo 12F IgG	Pneumococcal Serotype 12F IgG	85974-4
5197341	Pneumo 15A IgG	Pneumococcal Serotype 15A IgG	n/a
5197342	Pneumo 15C IgG	Pneumococcal Serotype 15C IgG	n/a
5197343	Pneumo 16F IgG	Pneumococcal Serotype 16F IgG	n/a
5575677	Pneumo 17F IgG	Pneumococcal Serotype 17F IgG	86010-6
5575703	Pneumo 19A IgG	Pneumococcal Serotype 19A IgG	40915-1
5197344	Pneumo 20A IgG	Pneumococcal Serotype 20A IgG	n/a
5575683	Pneumo 22F IgG	Pneumococcal Serotype 22F IgG	86053-6
5197345	Pneumo 23A IgG	Pneumococcal Serotype 23A IgG	n/a
5197346	Pneumo 23B IgG	Pneumococcal Serotype 23B IgG	n/a
5197347	Pneumo 24F IgG	Pneumococcal Serotype 24F IgG	n/a
5197348	Pneumo 31 IgG	Pneumococcal Serotype 31 IgG	n/a
5575709	Pneumo 33F IgG	Pneumococcal Serotype 33F IgG	47349-6
5197349	Pneumo 35B IgG	Pneumococcal Serotype 35B IgG	

<b>CPT Code(s)</b>	86581
<b>Lab Section</b>	Reference Lab