

## Eosinophil Count, Absolute

Order Name: **EOS CT ABS**

Test Number: 0100050

Revision Date: 10/31/2019

TEST NAME	METHODOLOGY	LOINC CODE
Eosinophil Count, Absolute	Flow cytometry	711-2

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	5 mL (1 mL)	Whole Blood	EDTA (Lavender Top)	Refrigerated
Alternate 1	1 mL (0.5 mL)	Whole Blood	EDTA (Lavender) Microtainer/Bullet	Refrigerated
Instructions	<b>For best results</b> Room temperature specimens should be tested within 24hrs, otherwise send Refrigerated. Refrigerated specimens can be tested up to 48hrs. Specimens received greater than 48hrs old will be canceled.			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1 Day
Clinical Use	Eosinophilia is found in parasitic diseases, allergic diseases, skin disorders, and certain malignancies.
CPT Code(s)	85004
Lab Section	Hematology