Specific Gravity Fluid

Order Name: SP GRAV FL Test Number: 0800900 Revision Date: 10/31/2019

TEST NAME		МЕТН	ODOLOGY	
Specific Gravity Fluid		Optica	d	2964-5
SPECIMEN REQUIRE	MENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	3mL (1 mL)	Body Fluid	EDTA (Lavender Top)	Ambient / Refrigerated
Instructions	Note fluid type on Tube and Request. Good for all Body Fluids except for CSF, Semen and Synovial. Invert EDTA tube several times to mix well and deliver to laboratory immediately. Specimen stability is 24hrs refrigerated. Rejection criteria: Specimens received in Syringes, urine cups, monovette tubes, red top tubes, and gel-barrier tubes. Frozen or clotted specimens will be rejected.			

GENERAL INFORMATION		
Testing Schedule	Daily	
Expected TAT	1 Day	
CPT Code(s)	84315	
Lab Section	Hematology	