

Adenosine Deaminase, Pericardial Fluid

Order Name: **Adeno D Pericard**

Test Number: 6904899

Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Adenosine Deaminase, Pericardial Fluid	Quantitative Spectrophotometry	49760-2

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.5 mL (0.2 mL)	Pericardial Fluid	Sterile Screwtop Container	Frozen
Instructions	Collect Pericardial Fluid, centrifuge specimen at room temperature. Transfer 0.5mL(0.1) fluid to a Standard Transport Tube and freeze. Indicate source on requisition. Specimen must remain frozen until received in lab. Unacceptable Conditions: Whole blood or Bronchoalveolar Lavage (BAL) specimens. Turbid specimens. Special Instructions: Testing referred to ARUP TC 3002984 Stability: Ambient 24 hours; Refrigerated 7 days; Frozen 30 days.			

GENERAL INFORMATION	
Expected TAT	8 - 13 days
Performing Labcorp Test Code	830932
Notes	Labcorp Test Code: 830932
CPT Code(s)	84311
Lab Section	Reference Lab