

Prenatal Profile #8 with HIV and HCV

Order Name: **Prenatal Pr8HC**
 Test Number: 6906279
 Revision Date: 03/19/2020

TEST NAME	METHODOLOGY	LOINC CODE
Complete Blood Count (CBC) with Automated Differential	Flow cytometry	See Individual Assays
Hepatitis B Surface Antigen	Chemiluminescence Assay	5195-3
Syphilis Antibody Screen	Chemiluminescence Assay	24110-9
Rubella Virus Antibody IgG		
HIV Antigen/Antibody Screen, 4th Gen	Immunoassay (IA)	56888-1
ABO Group & Rh Type		
Antibody Screen (Indirect Coombs)		
Hepatitis C Antibody (HCV Ab)	Chemiluminescence Assay	16128-1

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	See Instructions	See Instructions	See Instructions	
Instructions	This profile requires the collection of several specimen types. Please collect one of each of the following specimens: One 7mL EDTA Pink top. One 5mL EDTA Lavender top. One 5mL Clot Activator SST tube.			

GENERAL INFORMATION	
Testing Schedule	Test dependent
Expected TAT	2 - 4 days
CPT Code(s)	85025, 87340, 86780, 86762, 86900, 86901, 86850
Lab Section	See Individual Assays