Propoxyphene, MS, Ur, MAT

Order Name: Propoxyphene UR

Test Number: 5195585 Revision Date: 01/09/2024

| IESI NAME | | | METHODOLOGY | LOINC CODE |
|---------------------------|--|------------------------|---|-----------------------|
| Propoxyphene, MS, Ur, MAT | | Mass Spectrometry (MS) | | |
| SPECIMEN REQU | IREMENTS | | | |
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 30 mL (15 mL) | Urine, Random | Sterile Urine Container Preservative Free | Room Temperature |
| Instructions | Specimen: 30mL(15mL) Random Urine Preservative free Urine Transport tube Room Temperature or Refrigerated. Cause for Rejection: Insufficient volume, no ID on container, urine from preservative tube. SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. | | | |

| GENERAL INFORMATION | | | |
|---------------------------------|---------------|--|--|
| Expected TAT | 5-8 days | | |
| Performing Labcorp Test Code | 701746 | | |
| CPT Code(s) | 80367 | | |
| Lab Section | Reference Lab | | |

Stability: Room Temperature 5 days, Refrigerated 14 days, Frozen 14 days.

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025